




Improvement of Tomo-PIV Analyses in the Nasal Cavities

S. M. Tauwald^{1,2}, F. Erzinger², M. Quadrio³, M. Rütten⁴, C. Stemmer⁵, and L. Krenkel^{1,2}

- 1) Regensburg Center of Biomedical Engineering, OTH and University of Regensburg, and Regensburg Center of Health Sciences and Technology, Regensburg, Germany
- 2) Technical University of Applied Sciences (OTH) Regensburg, Department of Biofluid Mechanics, Regensburg, Germany
- 3) Politecnico Milano, Department Aerospace Science and Technology, Milan, Italy
- 4) German Aerospace Center (DLR), Institute of Aerodynamics and Flow Technology, Göttingen, Germany
- 5) Technical University of Munich (TUM), Department of Aerodynamics and Fluid Mechanics, Munich, Germany



Overall objective: bringing CFD into the clinical setting

↓
Patient-specific CT-scan

CFD solution

Implementation into surgical treatment


CFD enables the improvement of:

- Diagnostic precision
- Successful surgical outcomes




Checking the reliability of the numerical model through a solid benchmark by **Experimental Investigations**

Lateral nasal wall, sagittal section, medial view [1]




Paranasal sinuses, coronal section, anterior view [1]

Section A - A






Highly complex, inter-individually varying geometry

[1] Gilroy, A. M., MacPherson, B. R., Ross, L. M., Broman, J., & Josephson, A. (Eds.). (2008). *Atlas of anatomy* (pp. 356-450). Stuttgart: Thieme.



- 800 L **fish tank** containing three portholes
- **3-axis traversing unit** (isel Germany AG)
- **CCD cameras** (Imager Pro X 2A, 1600x1200 pix)
 - Macro lenses (Milvius 2/135, ZEISS)
 - Scheimpflug adapters
 - Cut-off filters
- **Light source** dp Nd:YAG laser (EverGreen200, 70-200 mJ @ 532 nm, 15 Hz)
- Trigger signal from control unit (PLC) starts **linear-motor-driven piston pump** → induces modified **breathing cycle**
- **Fluorescence tracer particles** (PMMA-RHB, 20-50 μm , $\text{abs/emm} = 560/584 \text{ nm}$)
- Two **seeding pumps**
- **Refractive Index Monitoring** laser and camera





Non-free-of-sight tomo-PIV requires Refractive Index Matching:




- $n_{(\text{water})} = 1.333$
 - $n_{(\text{glycerol})} = 1.473$
 - $n_{(\text{PDMS})} = 1.414$ (material data sheet)
- } working fluid

Problem:


- Customers value changed during processing
- Separation of the working fluid's components






Unsuccessful VSC results !



Laser visualization moves due to RI mismatch ing different RIs Refractive Index




$\dot{V}_1(\text{BCP})$ = Flow rate of seeding pump (Breathing Cycle Phase dependent)
 $\Delta t_2(\text{BCP})$ = Seeding time of seeding pump (Breathing Cycle Phase dependent)




FoV = 49.5 mm x 36.5 mm
RoI = 49.5 mm x 36.5 mm x 4.5 mm
Scale factor = 29.1 pix/mm
VSC error < 0.1 pix



Results – Averaging Cycle Phases




Averaged 3D Velocity Field





Velocity IVI / m/s






Snapshot of standard deviation of exemplary slice (IVI / m/s)







Velocity $|\mathbf{V}|$ / m/s





PIV-dt = 7 ms

Interrogation window: 32 vox








Frontal view





20.7 sec

Side view



21.6 sec




Successful experimental parts:

- **Phantom model development** method enables to build a patient-specific measurement geometry
- **In-site RIM method** minimizes errors and leads to acceptable Volume Self Calibration results
- **Linear-motor-driven-piston pump** gives the possibility to create **physiological flow rates**
- Analysing **different flow rates during cycle** with low-speed setup works while adapting **PIV-dt**
- **Triggered seeding strategy** leads to successful **seeding densities**
- **Tomo-PIV** enables investigations of **complex 3D flow structures**

Necessary improvements and further investigations:

- **High-speed system** will allow for a breathing cycle phase dependent PIV-dt
- **Seeding strategy optimization** for low flow regimes (Sinuses)
- Avoid separation of **working fluid's components**
- Enlargement **Field of View** -> Current: Head requires more than 200 RoIs
- Minimize **phase-locked shift** of 120 us



Thank you !



Contact

Sandra Melina Tauwald
OTH Regensburg
Dept. of Biofluid Mechanics
melina.tauwald@oth-r.de

Acknowledgments

This work is part of the cooperation with Politecnico di Milano and DLR Germany as part of the OpenNOSE project group

Funding

This work is supported by the Bavarian Academic Forum (BayWISS) – Doctoral Consortium “Health Research”, the Regensburg Center of Biomedical Engineering (RCBE), and Regensburg Center of Health Sciences and Technology (RCHST).